



# Cannabis Potency Results

Date: 11-Apr-15, 13:50:57

Client:


Sample: 70 Micron

Type:

Notes:

www.TestedLabs.com 207-618-9333  
(http://www.TestedLabs.com)

Cannabinoid	Result
THC	1.37 %
THCA	45.85 %
CBD	0 %
CBDA	0.81 %
CBN	0.17 %
CBG	0 %
THCV	0.34 %
CBC	0 %



<b>Total CBD</b>	0.81 %
<b>Max CBD</b>	0.71 %
<b>Total THC</b>	47.22 %
<b>Max THC</b>	41.59 %
<b>Total Cannabinoid Results</b>	48.54 %

- TestedLabs has tested your sample on a calibrated high performance liquid chromatography (HPLC) system. This method of testing allows us to measure both the acid (raw) and neutral (decarboxylated) states of the major cannabinoids in medical cannabis flowers and products.
- Liquid sample results are reported in milligrams per milliliter (mg/ml). All other samples are reported as percent per dry weight.
- 1ml = approx 1 standard dropper full or 20 drops. 5 ml = 1 teaspoon.
- Please use a graduated dropper or oral syringe for accurate dosing.
- Please follow the dosage recommendations made by your medical provider. Max THC and Max CBD are estimations of the total amounts of each cannabinoid based on 100% decarboxylated THCa and CBDa. (Using the multiplier 314.468/358.471).
- Heat activation of cannabis products converts THCA to THC and CBDA to CBD in a time and temperature dependent manner. THCa typically converts fully to THC at 252 Fahrenheit for 27 minutes.

Reviewed by Clay Sulak \_\_\_\_\_ on 04/14/2015